

# Bead Ruptor 12

## BEAD MILL HOMOGENIZER



*Compelling performance  
at a compelling price*



# Bead Ruptor 12



ORDER NO. 19-050A

*Ideally suited for nucleic acid, drug, protein, enzyme, and virus extractions.*

## Specifications

Processing Range: 25  $\mu$ L to 6 mL

Power Rating: 115-230 VAC/50-60 Hz

Dimensions: 13" W x 15.5" D x 13" H

Weight: 29 kg (65 lbs.)

Compliance: CE

Run Time: 1s - 9:59 m

Performance Range: 0.8 m/s - 6 m/s  
in increments of 0.15 m/s

Cycles: 1-10


## Includes Tube Carriage for 2 mL and 7 mL tubes

- Holds 12 x 2 mL or 0.5 mL tubes and/or 4 x 7 mL tubes
- Includes semi-permanent ring



## Features

- Optimized tube motion with speeds up to 6 m/s
- Process up to 16 samples per cycle
- Convenient front loading design
- Locking safety lid
- Reproducible extraction of proteins, nucleic acids, and small molecules



935C Cobb Place Blvd.  
Kennesaw, GA 30144  
800.776.4431 • 770.421.0058  
www.omni-inc.com



REV.A