

Thermocycler	Thermal-Cycler	QPCR-Cycler	Block								
				s	s	s	lp	lp	s	s	
standard (s) oder low profile (lp)											
8er-Streifen/FrameStrip® 8er-Strips/Strips-Caps/PCR-Tubes				34-210-0125	34-220-0125	34-224-0120	34-225-0120	34-228-012	34-230-xxx	34-260-xxx	
ABI/Life Technologies/Thermo Fisher											
Veriti, Proflex, Simplicamp	•		96 well	✓	✓	✓			✓	✓	
GeneAmp® 2700/2720/9600/9700	•		96 well	✓	✓	✓			✓	✓	
GeneAmp® 9800 Fast, Veriti FAST	•		96 well								
GeneAmp® 9700, Veriti, Proflex, Multiblock System	•		96 well								
7000,7300,07500,7700,7900HT		•	96 well	✓	✓	✓			✓	✓	
QuantStudio 3/5/6/7/12K, ViiA7		•	96 well	✓	✓	✓					
StepOne		•	96 well					✓			
StepOnePlus, 750FAST, 7900HTFAST, QuantStudio, Fast		•	96 well				✓	✓			
QuantStudio 3/5/6/7/12K, ViiA7 Fast			96 well	✓	✓	✓	✓	✓			
Sequencers 3100	•		96 well								
Sequencers	•		96 well								
Agilent/Stratagene											
Surecycler 8800	•		96 well	✓	✓						
Robocycler Gradient	•		96 well				✓	✓			
AriaMx		•	96 well				✓	✓			
Mx3000P		•	96 well	✓	✓		✓	✓	✓	✓	
Mx4000P		•	96 well	✓	✓				✓	✓	
Analytic Jena/Biometra											
TRIO, Tpersonal, T3 Thermocycler	•		96 well	✓	✓	✓			✓	✓	
Flexcycler2, T1, Tgradient, Tone, Tadvanced, Tprofessional	•		96 well	✓	✓	✓			✓	✓	
Trobot 96/SpeedCycler2	•		96 well	✓	✓	✓	✓	✓	✓	✓	
qTOWER/G/touch/Toptical		•	96 well	✓	✓	✓			✓	✓	
BIOER Technologies											
GeneQ	•		96 well	✓	✓	✓			✓	✓	
Gene Touch 96	•		96 well	✓	✓	✓			✓	✓	
Bio-Rad											
Genecycler	•		96 well	✓	✓	✓			✓	✓	
C1000 Touch, S1000	•		96 well	✓	✓	✓	✓	✓			
iCycler™, MyCycler™, T100	•		96 well	✓	✓	✓			✓	✓	
CFX96Touch, CFX96DeepWell Connect		•	96 well				✓	✓			
MyiQ, iCycler IQ		•	96 well				✓	✓	✓	✓	
Bio-Rad/MJ Research											
Mini Gradient	•		96 well	✓	✓	✓			✓	✓	
Personal	•		96 well								
PTC100™, 200™, 220™, 221™, 225™, 240™	•		96 well	✓	✓	✓	✓	✓	✓	✓	
MiniOpticon		•	96 well				✓	✓			
Chromo4		•	96 well	✓	✓	✓	✓	✓	✓	✓	
Opticon2		•	96 well				✓	✓			
Corbett Research											
(Qiagen) Palm Cycler	•		96 well	✓	✓	✓	✓		✓	✓	
qPCR Cycler		•	96 well								
Eppendorf											
Mastercycler® ep, Pro, Pro S, nexus	•		96 well	✓	✓	✓	✓	✓	✓	✓	
Mastercycler® nexus X2, GX2, GX2e, X2e	•		96 well	✓	✓	✓	✓	✓	✓	✓	
Mastercycler® ep realplex		•	96 well	✓	✓	✓	✓	✓	✓	✓	
GE Healthcare/Amersham											
Sequencer: MegaBACE™500, MagaBACE™1000 mark2			96 well					✓			
Sequencer: MegaBACE™4000			96 well					✓			
PEQLAB/VWR											
peqSTAR XS, peqSTAR 2X	•		96 well	✓	✓	✓	✓	✓	✓	✓	
peqSTAR 96X	•		96 well	✓	✓	✓	✓	✓	✓	✓	
ROCHE		•									
LC96, LC480			96 well					✓			
Sensoquest											
Labcycler	•		96 well	✓	✓	✓	✓	✓	✓	✓	
Takara											
Dice touch, Gradient	•		96 well	✓	✓	✓			✓	✓	
Techne											
3Prime, 3PrimeG	•		96 well	✓	✓	✓			✓	✓	
Prime, PrimeG, Prime Elite, Prime Elite Satellite	•		96 well	✓	✓	✓			✓	✓	
PCRmax Alpha Cycler 1/2/4	•		96 well	✓	✓	✓			✓	✓	
TC412, TC512, Genius, Genius Quad, Touchgene, Flexigene	•		96 well	✓	✓	✓	✓	✓	✓	✓	
Quanta		•	96 well				✓	✓			

✓ kompatibel

Stand: Januar 2023